

LUSTER[®]-Kits for Zirconia

MEISINGER CAD/CAM polishing system for zirconium oxide

The instruments of this LUSTER[®] Set were developed specifically for working with zirconium oxide ceramics. By working the material with these special polishers, the crown attains its maximal gloss.



LUS85

Adjusting & Polishing Kit (LUS85)

Kit with grinding and polishing instruments (conical form, wheel and flame) for CAD/CAM restorations for specific treatment of monolithic zirconium oxide, for example, BruxZir[®] crowns. The instruments are also suited as correction tools for glass corrosion defects. Convincing outcome with exceptional gloss. All instruments for use with the handpiece.



LUS91

Intraoral Adjustment Kit (LUS91)


The kit consists of three specifically developed diamonds for zirconium oxide and two each of perfectly matched polishers interspersed with diamonds (cup and point-shaped flame form) for smoothing and pre-polishing and for high-gloss polishing of zirconium oxide ceramics. The instruments are for intraoral use following placement of crowns in patients. Instruments with FG or RA shank.

Our sales force would be pleased to answer your questions: +49 2131 2012-303

Meisinger
125
YEARS

LUSTER®-Kits for Zirconia

Fax for ordering +49 2131 2012-222

	Figure		Size	Grain	recomm. / max. revolutions	Packaged à	Price / Piece*	Amount 
	LUS85	LUSTER® for Zirconia	-	-	-	1	159,- €	
	9736-HP	Pre-Grinding	150	coarse	8.000 - 12.000 / 25.000 rpm	2	29,90 €	
	Z652R-HP	Grinding	035	medium	10.000 - 15.000 / 30.000 rpm	2	15,90 €	
	DCA04-HP	Smoothing and pre-polishing	040	medium	7.000 - 12.000 / 20.000 rpm	2	10,90 €	
	DCA06	Smoothing and pre-polishing	170	medium	7.000 - 12.000 / 20.000 rpm	1	29,90 €	
	DCA10-HP	High gloss polish	040	fine	7.000 - 12.000 / 20.000 rpm	2	10,90 €	
	DCA12	High gloss polish	170	fine	7.000 - 12.000 / 20.000 rpm	1	29,90 €	
	LUS91	LUSTER® for Zirconia intraoral	-	-	-	1	79,- €	
	Z801-FG	Diamond for zirconium oxide	010	medium	160.000 / 300.000 - 350.000 rpm	5	5,40 €	
	Z850-FG	Diamond for zirconium oxide	018	medium	160.000 / 300.000 - 350.000 rpm	5	5,40 €	
	Z833-FG	Diamond for zirconium oxide	023	medium	160.000 / 300.000 - 350.000 rpm	5	5,40 €	
	DCA04-RA	Smoothing and pre-polishing	040	medium	7.000 - 12.000 / 20.000 rpm	2	10,90 €	
	DCA13-RA	Smoothing and pre-polishing	060	medium	7.000 - 12.000 / 20.000 rpm	2	10,90 €	
	DCA10-RA	High gloss polish	040	fine	7.000 - 12.000 / 20.000 rpm	2	10,90 €	
	DCA14-RA	High gloss polish	060	fine	7.000 - 12.000 / 20.000 rpm	2	10,90 €	

* plus VAT

CE 0044

Available from the dental trade.

Your trade depot:

Your delivery address:

Hager & Meisinger GmbH

Hansemannstr. 10 • 41468 Neuss • Germany
 Phone: +49 2131 2012-0 • Fax: +49 2131 2012-222
 info@meisinger.de • www.meisinger.de

Meisinger USA, L.L.C.

10200 E. Easter Avenue • Centennial, Colorado 80112 • USA
 Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
 Toll free: +1 (866) 634-7464
 info@meisingerusa.com • www.meisingerusa.com

85FLO99E-0213